

IMAGING ELEVATED


9800 Efficiency Kit



Includes:

- **Wireless footswitch control**
- **High Definition LED Touchscreen Display**
- **Wireless DICOM + secure USB to GE OEC 9800 C-Arms**

Make your 9800 images as good as a 9900 for 1/10th of the price!



OEC Monitor Upgrade replaces CRT's with touch-screen LED

Flash Adds wireless DICOM + secure USB



Spectre

Wireless Footswitch



Spectre Wireless Footswitch

Cut the cable, get the freedom!

- Compatible with GE-OEC Series 9800 C-Arms
- Plug and play connectivity
- Encrypted transmission
- Factory-mated units - Invisible to other devices
- FDA 510k & ETL approved
- Will not interfere with hospital equipment.
- Range: up to 50 feet



OEC Monitor Upgrade

- Touchscreen high resolution LED display
- New monitor bezel
- Mounting plates
- Step-by-step installation instructions

**Approved for all medical safety standards
FDA registered**

Two year full replacement warranty



Flash Wireless DICOM and Secure USB

Lose the hassle, gain connectivity!

- Adds up to 10,000 new image storage
- Secure access controlled USB
- Connectable to any system with DICOM
- Fully compatible with the LED monitor upgrade (installed as part of the Monitor upgrade).
- Includes all available enterprise security software

