

# E-CUBE 15 EX

Product Specification

V3.2 Rev.01



**ALPINION**  
MEDICAL SYSTEMS

## General Specifications

### Physical Specification

- Width: 580 mm
- Depth: 895 mm
- Height: 1,413 - 1,848 mm
- Weight: 105 Kg (System only)

### System Design

- 19.5" LCD monitor and articulated arm
- Front and Rear handle
- 4 active probe connectors
- 1 pencil probe port
- 4 swivel wheels
- Integrated HDD
- Built-in DVD-R drive
- ECG (Optional)
- Thermal printer storage
- Built-in stereo speaker
- Built-in cable management
- Anti-dust probe connector door
- System boot-up: 120sec

### Main Monitor

- 19.5" high definition LED
- Resolution: 1366 X 768 (16:9 wide)
- 16.2M color display
- Contrast: 1:1,000
- **IPS (In Plane Switching)** technology
- Bright and contrast adjustment
- Articulated monitor arm
  - Up/down: 140 mm
  - Swivel: +/- 90°
  - Tilt: +15° / -90°

### Control Panel

- **Electric** motorized control panel adjustment
  - Up/down: 175 mm
  - Swivel: +/- 25°
- Alphanumeric keyboard with backlit
- Ergonomic key layout
- 8 TGC control levers
- 6 user customizable keys
- 5 removable probe holders
- **Built-in gel warmer**

### Touch Screen

- 10.4" high definition LED
- Resolution: 1024 X 768
- 16.2M color display
- Bright and contrast adjustable (OSD)

### Computing Power

- Main processor: Intel core i5-2510E 2.5 GHz
- Main memory: 4 GB
- HDD: 500 GB

### Electrical Specification

- 100 - 120 V, 220-240 V, 50/60 Hz
- Power consumption: Max. 900 VA

## System Standard Feature

### Architecture

- Digital beamforming
- High density transducer connector
- **Triple** multi-beam
- Flexcan™
- SensitiView™
- Crystal Signature™
- MicroFit™
- Optimal Imaging Suite™
- Needle enhancement
- **Multi-layer transducer**
- OS: WindowsXP Embedded

### Transducer Type

- Convex
- Micro convex
- Endo cavity
- Linear
- Phased (Sector)
- Volume
- Pencil
- TEE

### Application

- Abdomen
- Breast
- Cardiac
- EM
- Gynecology
- Musculoskeletal
- Obstetrics
- Pediatric
- Small Parts
- Urology
- Vascular

### Preset

- Abdomen
- Abdomen (Pen)
- Renal
- OB1, OB2, OB3

- Fetal Echo
- Thyroid
- Testicle
- MSK
- Superficial
- MSK deep
- Nerve
- Shoulder/Knee
- Wrist/Elbow
- Breast
- EM
- Cardiac
- Cardiac (Pen)
- TCD

### Imaging Mode

- 2D (B)
- Color flow (CF)
- Pulse wave doppler (PW)
- Power doppler (PD)
- Directional power doppler (DPDI)
- Tissue doppler Imaging (TDI)
- Continuous wave Doppler (CW)
- Motion (M)
- Color M
- Anatomical M
- 3D / 4D (Volume)
- Live HQ™
- Panoramic
- Elastography
- Needle Vision™ / Needle Vision™ Plus
- Dual Live

### Complex Mode

- Duplex: B+B, B+PW, B+CW, B+CF, B+PD, B+DPDI, B+M, B+3D, B+4D, B+Live HQ, A+B+C+3D, A+B+C+4D, B+Elastography
- Triplex: B+CF+PWD, B+CF+CWD, B+CF+M

### Scanning Parameters

- System frequency range: 1 ~ 20 MHz
- Processing channel: 294,912
- ADC: 12 bits
- Displayed imaging depth: 0 ~ 30cm
- Max. frame rate (Probe dependant)
  - 2D: 1,369 (Hz/FPS)
  - Color: 269 (Hz/FPS)
  - Volume: 70 (Hz/VPS)
- Dynamic range: Max. 192dB
- Gray scale: 256 levels
- Focus
  - Focal number: Max.8
  - Focal position: Max. 17

- Zoom
  - Write zoom
  - Read zoom
- CF
  - 300 Hz ~ 18.1 KHz
  - 6 cm/s ~ 34.8 m/s
- PW
  - 330 Hz ~ 18.1 KHz
  - 2 cm/s ~ 50.4 m/s
- CW
  - 300 Hz ~ 78.1 KHz
  - 2.5 cm/s ~ 300.7 m/s

### Imaging Presentation

- Dual
- Quad
- Duplex image format
  - 1:1 / 1:2 / 2:1 / full screen
- Image reverse: Left/right
- Image rotation: 0°, 90°, 180°, 270°

## Display and User Interface

### Display Annotation

- Patient ID: Up to 64 characters
- Patient name: First, last and middle
- Gender, age and birth date
- Gestational age from LMP, EDD and GA
- Transducer name
- Institution / hospital name: up to 25 characters
- Date format: 3 types
  - YYYY/MM/DD
  - MM/DD/YYYY
  - DD/MM/YYYY
- Time format: 2 types
  - 12 hours
  - 24 hours
- Acoustic power
  - Mechanical index (MI)
  - Thermal index (TI): TIS, TIB, TIC
- Transducer name
- Preset
- Transducer orientation
- Depth
- Focus
- TGC line
- Frequency
- Frame rate
- Power
- Mode name
- 2D mode
  - Gain
  - Dynamic range

- Gray scale bar
- Doppler mode
  - Wall filter
  - Persistence
  - PRF
  - Sample volume position, size and angle
  - Velocity or frequency scale
  - Time
- Color flow mode
  - Wall filter
  - ROI box
  - PRF
  - Velocity or frequency scale
  - Color scale bar
- M mode
  - Gain
  - Dynamic range
  - Time
- Arrow
- Indicator
- Text
  - Size: Small, medium, large
- Body pattern (Marker)
  - 201 sets
  - Size: Small, medium, large
- Biopsy guide line
- Thumbnail
- Stored image info
- Cine gauge
- Measurement result

## Live Mode Display

- 2D
  - Angle steer (Probe dependent): 7 steps (-15 ~ 15°)
  - Colorize: 0 ~ 254
  - Dynamic range: 30 ~ 192dB
  - Edge enhance: 0 ~ 4
  - Frequency: 3 steps
  - Full SRI: 0 ~ 5
  - Gray map: 0 ~ 37
  - Gain: 0 ~ 100
  - Harmonic
  - Pulse inversion harmonic
  - Invert
  - Line density: 0 ~ 5
  - Persistence: 0 ~ 5
  - Power: 1 ~ 100%
  - Reject: 1 ~ 9
  - Spatial compound: 0 ~ 3
- M
  - Anatomical M: On, off
  - AMM angle
  - Colorize: 0 ~ 22

- Dynamic range: 30 ~ 150dB
- Gray map: 0 ~ 22
- Gain: 0 ~ 100
- Reject: 0 ~ 10
- Power: 1 ~ 100%
- Sweep speed: 0 ~ 6
- PWD
  - Angle: -89° ~ 89°
  - Angle steer (Probe dependent): 7 steps (-15 ~ 15°)
  - Auto angle
  - Base line: -8 ~ 8
  - Dynamic range: 30 ~ 120dB
  - Frequency: 3 steps
  - Gray map: 0 ~ 13
  - Gain: 0 ~ 100
  - Line density: 0 ~ 5
  - Power: 1 ~ 100%
  - Reject: 0 ~ 10
  - Scale (PRF): 0.3 ~ 18.1 KHz
  - SV size: 0.7 ~ 15mm
  - Sweep speed: 0 ~ 6
  - Sound: 0 ~ 100%
  - Wall filter: 0 ~ 7
- CWD
  - Angle: -89° ~ 89°
  - Base line: -8 ~ 8
  - Dynamic range: 30 ~ 120dB
  - Gray map: 0 ~ 13
  - Gain: 0 ~ 100
  - Scale (PRF): 0.3 ~ 78.1 KHz
  - Frequency: 3 steps
  - Power: 1 ~ 100%
  - Reject: 0 ~ 10
  - Sound: 0 ~ 100%
  - Sweep speed: 0 ~ 6
  - Wall filter: 0 ~ 9
- CF/PD
  - Angle steer (Probe dependent): 7 steps (-15 ~ 15°)
  - Auto angle
  - Base line: -8 ~ 8
  - Color map: 13
  - Dynamic range: 30 ~ 120dB
  - Ensemble: 6 ~ 16
  - Frequency: 3 steps
  - Flow state: Low, mid, high
  - Gain: 0 ~ 100
  - Line density: 0 ~ 5
  - Power: 1 ~ 100%
  - Scale (PRF): 0.3 ~ 18.1 KHz
  - Smooth: 0 ~ 10
  - Wall filter: 0 ~ 7
- 3D/4D
  - Quality: Low, Mid, H1, H2

- Volume angle: 15 ~ 75°
- View direction: Up, down, Lt, Rt, front, back
- Threshold: 0 ~ 255
- Brightness: 1 ~ 100
- Contrast: 1 ~ 100
- MPR
  - Display format
  - Ref. plane
  - 3D orientation
  - Smooth: 1 ~ 6
  - Render mode: Surface, MaxIP, MinIP, XRay, Light
  - Volume review
  - Easy cut
- Cube CT
  - Play / Stop
  - Play mode: YoYo, Loop
  - Rotation Axis: X, Y, Z
- Multi Slice
  - Display format: 1\*1 / 1\*2 / 2\*2 / 3\*3 / 4\*4
  - Reference Plane: A / B / C
- Live HQ
  - Light Direction
  - HQ Scatter: 0 ~ 10
  - Move light
- Panoramic
  - Colorize: 0 ~ 254
  - Dynamic range: 30 ~ 192dB
  - Edge enhancement: 0 ~ 4
  - Frequency: 3 steps
  - Full SRI: 0 ~ 5
  - Gray map: 0 ~ 37
  - Gain: 0 ~ 100
  - Harmonic
  - Pulse inversion harmonics
  - Line density: 0 ~ 5
  - Persistence: 0 ~ 5
  - Power: 1 ~ 100%
  - Reject: 1 ~ 9
- Elastography
  - Live dual: On / off
  - Invert: On / off
  - E persistence: 1 ~ 4
  - Alpha blending: 0 ~ 100
  - Color map: 0 ~ 8
- Needle Vision plus
  - L/R Flip
  - Angle: Shallow, Medium, Steep

## Post Processing

- 2D
  - Dynamic range
  - Colorize
  - Full SRI / SRI

- Gray map
- Reject
- TGC
- M
  - AutoCalc
  - Baseline
  - Colorize
  - Dynamic range
  - Gray map
  - Reject
  - Sweep speed
- PWD
  - Angle correction
  - AutoCalc
  - Baseline
  - Colorize
  - Dynamic range
  - Gray map
  - Reject
  - Sweep speed
- CWD
  - Angle correction
  - AutoCalc
  - Baseline
  - Colorize
  - Dynamic range
  - Gray map
  - Reject
  - Sweep speed
- CF/PD
  - Baseline
  - Colorize
  - Color map
  - Threshold

## Technology and Feature

### Imaging Technology

- Filtered tissue harmonic imaging (FTHI)
- Pulse inversion tissue harmonic imaging (PTHI)
- Speckle reduction imaging (SRI)
- Full speckle reduction imaging (Full SRI)
- Frequency compounding imaging (FCI)
- Needle Vision™
- Needle Vision™ Plus
- Xspeed™
- Live HQ™
- Beam steer
- Virtual Convex (Trapezoidal)
- S-FOV (Extended view)
- Automated intima media thickness measurement

## Features (Options)

- Auto NT
- Auto IMT Measurement
- CV measure & report
- SRI
- Full SRI
- Xpeed™
- Frequency Compounding
- Spatial Compounding
- Filtered Tissue Harmonic Imaging (FTHI)
- Continuous Wave (CW)
- DICOM
- DICOM SR: OB, Adult Echo, Vascular
- Anatomical M
- Panoramic
- Elastography
- Tissue Doppler Imaging (TDI)
- Live HQ™
- Cube View™
- Cube Strain™
- Stress Echo
- Needle Vision™
- Needle Vision™ Plus
- 3D/4D
- 4 probe port
- ECG (IHE / AHA)
- Pencil probe port
- Side basket
- Foot switch (3 pedals)

## Transducer

### SC1-4H

- Convex array
- Frequency: 1.0 - 4.0 MHz
- Field of view: 60°
- Elements: 192
- Radius of curvature: 60 mm
- Footprint: 72.4 X 16.8 mm
- Application: Abdomen, EM, Gynecology, Obstetrics
- Crystal Signature Technology
- Biopsy available

### SC1-4HS

- Convex array
- Frequency: 1.0 - 4.0 MHz
- Field of view: 73°
- Elements: 166
- Radius of curvature: 46 mm
- Footprint: 68 X 15.9 mm
- Application: Abdomen, EM, Gynecology, Obstetrics
- Crystal signature technology

- Biopsy is not available

### SC1-6H

- Convex array
- Frequency: 1.0 - 6.0 MHz
- Field of view: 60°
- Elements: 192
- Radius of curvature: 60 mm
- Footprint: 72.4 X 16.8 mm
- Application: Abdomen, EM, Gynecology, Obstetrics
- Crystal signature technology
- Biopsy available

### C5-8

- Micro convex array
- Frequency: 5.0 - 8.0 MHz
- Field of view: 92°
- Elements: 128
- Radius of curvature: 14 mm
- Footprint: 35 X 6.2 mm
- Application: Abdomen, Cardiac, EM
- Biopsy is not available

### C5-8N

- Micro convex array
- Frequency: 5.0 - 8.0 MHz
- Field of view: 93.6°
- Elements: 128
- Radius of curvature: 15 mm
- Footprint: 25.5 X 9 mm
- Application: Abdomen, Cardiac, EM
- Biopsy is not available

### L3-8

- Linear array
- Frequency: 3.0 - 8.0 MHz
- Field of view: 38.4 mm
- Elements: 128
- Footprint: 44.8 X 7.8 mm
- Application: Breast, EM, MSK, Vascular, Small Parts
- Biopsy available

### L3-12H

- Linear array
- Frequency: 3.0 - 12.0 MHz
- Field of view: 38.4 mm
- Elements: 192
- Footprint: 44.8 X 7.8 mm
- Application: Breast, EM, MSK, Vascular, Small Parts
- Biopsy available

### L3-12X

- Linear array
- Frequency: 3.0 - 12.0 MHz
- Field of view: 51.2 mm
- Elements: 256
- Footprint: 58.2 X 10.4 mm
- Application: Breast, EM, MSK, Vascular, Small Parts
- Biopsy is not available

### L8-17X

- Linear array
- Frequency: 8.0 - 17.0 MHz
- Field of view: 51.2 mm
- Elements: 256
- Footprint: 58.2 X 10.4 mm
- Application: Breast, EM, MSK, Vascular, Small Parts
- Biopsy is not available

### IO3-12

- Linear array
- Frequency: 3.0 - 12.0 MHz
- Field of view: 16 mm
- Elements: 80
- Footprint: 22.8 X 6 mm
- Application: EM, Small Parts
- Biopsy is not available

### IO8-17

- Linear array
- Frequency: 8.0 - 17.0 MHz
- Field of view: 25.6 mm
- Elements: 128
- Footprint: 36.7 X 4.2 mm
- Application: Small Parts, MSK
- Biopsy is not available

### SP1-5X

- Phase array
- Frequency: 1.0 - 5.0 MHz
- Field of view: 90°
- Elements: 96
- Footprint: 24.8 X 17.6 mm
- Application: Abdomen, Cardiac, EM, TCD
- Crystal signature technology
- Biopsy is not available

### MP1-5X

- Phased array
- Frequency: 1.0 - 5.0 MHz
- Field of view: 90°

- Elements: 96
- Footprint: 24.8 X 17.6 mm
- Application: Abdomen, Cardiac, EM, TCD
- Biopsy is not available

### SP3-8

- Phase array
- Frequency: 3.0 - 8.0 MHz
- Field of view: 90°
- Elements: 64
- Footprint: 16 X 12.8 mm
- Application: Abdomen, Cardiac, EM, Pediatric
- Crystal signature technology
- Biopsy is not available

### E3-10H

- Endo cavity
- Frequency: 3.0 - 10.0 MHz
- Field of view: 150°
- Elements: 192
- Footprint: 21.5 X 18.2 mm
- Application: GYN, OB, Urology, EM
- Biopsy available

### EC3-10H/EV3-10H

- Endo cavity
- Frequency: 3.0 - 10.0 MHz
- Field of view: 156°
- Elements: 192
- Footprint: 21.5 X 18.6 mm
- Application: GYN, OB, Urology, EM
- Biopsy available

### EC3-10X/EV3-10X

- Endo cavity
- Frequency: 3.0 - 10.0 MHz
- Field of view: 230°
- Elements: 256
- Footprint: 21.5 X 18.6 mm
- Application: GYN, OB, Urology, EM
- Biopsy available

### SVC1-6H

- Volume convex array
- Frequency: 1.0 - 6.0 MHz
- Field of view: 79°
- Sweep angle: 75°
- Radius of curvature: 40 mm
- Elements: 192
- Footprint: 59.2 X 45.2 mm
- Application: Abdomen, OB, GYN, EM
- Crystal signature technology

- Biopsy is not available

### VE3-10H

- Volume endo cavity
- Frequency: 3.0 - 10.0 MHz
- Field of view: 160°
- Sweep angle: 105°
- Radius of curvature: 10 mm
- Elements: 192
- Footprint: 24.4 X 24.4 mm
- Application: GYN, OB, Urology
- Biopsy is not available

### TEE3-7

- Phased array
- Frequency: 3.0 - 7.0 MHz
- Field of view: 90°
- Footprint: 11 X 11 mm
- Array rotation: 0 - 180°
- Application: Cardiac, EM
- 4 ways control
- Biopsy is not available

### CW2.0

- Pencil
- Frequency: 2.0 MHz
- Footprint: 17.2 X 17.2 mm
- Application: Cardiac
- Biopsy is not available

### CW5.0

- Pencil
- Frequency: 5.0 MHz
- Footprint: 17.2 X 17.2 mm
- Application: Cardiac
- Biopsy is not available

## Data Management

### Image Archive

- E-View™
- Image store: Up to 52,000 images (382 GB)
- Image format
  - BMP, JPEG, DICOM, AVI, WMV
- Export / Backup / Restore
- Network storage
- Cube View™
- DVR

### Cine

- Cine memory: 500 MB, 2,000 frames, 35 volumes

- Prospective and retrospective
- Cine review, save and edit
- Run, stop, frame move
- Play speed: 50% 100%, 200%, 300%, 400%

## Measurement and Calculation

### Basic (Caliper)

- 2D
  - Distance
  - Area circumference
  - Volume
  - Angle
  - A/B ratio
  - Disk Volume
- M
  - Distance, Time, Slope,
  - %Stenosis
  - HR
  - A/B Ratio
  - Distance
- Doppler
  - Velocity
  - Time
  - PI
  - Trace
  - Acceleration
  - A/B Ratio
  - VTI
  - RI
  - Volume Flow

### Abdomen Measurement

- AO
- Renal
- %Stenosis
- A/B ratio

### Breast

- 2D
  - Lesion
- M
  - HR
- Doppler
  - Ves
  - HR

### Cardiology Measurement

- 2D
  - Ao/LA
  - LA Vol(Diam)
  - Aorta Diam



- Pul. Diam
- Vena Cava
- RV
- Simpson BP
- Simpson SP
- Modified Simpson
- Area Length
- Teichholz
- LV Mass T-E
- LV Mass A-L
- LA Vol (A-L)
- LA Vol (Simpson)
- RA Vol (A-L)
- RA Vol (Simpson)
- TV
- PISA
- MV
- AV
- PV
- Shunts
- AVDia

- M

- Teichholz
- Ao/LA
- MV
- RV
- Pul. Vein

- Doppler

- TV
- Pul. Vein
- TDI
- AV
- PV
- Shunts
- PISA
- MV

## Emergency Medicine

- 2D

- AO
- Renal L
- Bladder Vol.

- M

- %Stenosis
- A/B Ratio
- HR

- Doppler

- AO
- Renal A
- SMA
- IMA
- IVC
- Uterine A

- Ovarian A
- Ves
- HR

## Gynecology Measurement

- 2D

- Cervix
- Endo
- OV Vol.
- UT Vol.
- Bladder Vol.
- OV
- UT
- Follicle
- Bladder
- Early Gest.
- Uterine

- M

- HR

- Doppler

- AO
- Umbilical A
- Placenta
- Des.Aorta
- Uterine A
- Mid Cerebral A
- Ovarian A
- Ves
- HR

## Obstetrics Measurement

- 2D

- AC
- HC
- FL
- BPD
- CRL
- GS
- MSD
- OFD
- YS
- Humerus
- Radius
- Ulna
- Fibula
- NB
- Foot
- MAD
- FAD
- TC
- APTD
- APAD
- APTDxTTD

- CTAR
- NT
- AFI
- M
  - HR
  - Heart
  - Teichholz
- Doppler
  - AO
  - Umbilical A
  - Placenta
  - Des. Aorta
  - Uterine A
  - MCA
  - Ovarian A
  - Ves
  - DV
  - Valves
  - Ventricles
  - Arteries
  - Veins
  - HR

### MSK Measurement

- 2D
  - Volume
  - Hip (BA) / Hip (AB)
- M
  - HR
- Doppler
  - Ves
  - HR

### Pediatric

- 2D
  - AO
  - Renal L
  - Renal Vol.
  - Hip
- M
  - HR
- Doppler
  - AO
  - Renal A
  - CHA
  - IVC
  - Hepatic V
  - Mid HV
  - MPV
  - PV
  - Splenic A
  - Splenic V
  - SMA

- SMV
- IMA
- Iliac A
- Iliac V
- Ves
- HR

### Small Parts

- 2D
  - Thyroid Vol.
  - Testicle Vol.
- M
  - HR
- Doppler
  - Ves
  - HR

### Urology

- 2D
  - Renal L
  - Renal Vol.
  - Bladder Vol.
  - Prostate Vol.
- M
  - HR
- Doppler
  - Ves.
  - HR

### Vascular

- Carotid 2D
  - CCA
  - ICA
  - Bulb
  - ECA
  - Vert. A
- Renal 2D
  - Arcuate A Sup.
  - Arcuate A Med.
  - Arcuate A Inf
  - Renal A
  - Kidney
- Aortic lilac. 2D
  - Renal A
  - CFA
- Mesenteric 2D
  - Abd. Ao
  - Celiac A
  - SMA
  - IMA
  - Hepatic A
  - Portal Vein
  - Splenic A

- Upper Extremity 2D
  - Subclav. A
  - Brachial A
  - Radial A
  - Ulnar A
- Lower Extremity 2D
  - CFA
  - SFA
  - PFA
  - Pop. A
  - PTA
  - ATA
  - Pero. A
- Carotid M / Renal M / Aorta lilac M / Mesenteric M / Upper Extremity M / Lower Extremity M
  - HR
- Carotid Doppler
  - CCA
  - ICA
  - ECA
  - Bulb
  - Vert. A
  - HR
- Renal Doppler
  - Arcuate A Sup.
  - Arcuate A Med.
  - Arcuate A Inf
  - Renal A
  - Kidney
  - HR
- Aortic lilac. M
  - Renal A
  - CFA
  - HR
- Mesenteric Doppler
  - Abd. Ao
  - Celiac A
  - SMA
  - IMA
  - Hepatic A
  - Portal Vein
  - Splenic A
  - HR
- Upper Extremity Doppler
  - Subclav. A
  - Brachial A
  - Radial A
  - Ulnar A
  - HR
- Lower Extremity 2D
  - CFA
  - SFA
  - PFA

- Pop. A
- PTA
- ATA
- Pero. A
- HR
- TCD
  - ICA
  - MCA
  - ACA
  - PCA P1
  - PCA P2
  - AComA
  - PComA
  - Vert A
  - BA
  - Ves

## Connectivity

### Peripheral

- USB 2.0 (8 ports)
- Ethernet 10/100 Base-T
- Remote control
- HDMI video / audio out
- DVI out
- VGA (D-Sub) out
- S-VHS out
- Composite video out
- BNC out
- Audio in L/R
- Audio out L/R

### DICOM

- DICOM 3.0
- Image storage
- Multi-frame image storage
- Secondary image capture storage
- Structured report
- Grayscale print management
- Color print management
- Storage commitment push model
- Modality performed procedure step
- Worklist
- Verification
- Media storage: USB media, DVD, CD-R

### Thermal printer

- BW
  - Sony UP-D897MD
  - Sony UP-D898MD
  - Sony UP-X898MD
  - Mitsubishi P95DE

- Color
  - Sony UP-D25MD
  - Sony UP-25MD
- WEEE
- Rohs

## Safety and Compliance

- IEC/EN 60601-1:6.8.1:
  - Type of protection against electrical shock: Class I
  - Degree of protection against electrical shock (Patient connection): Type BF equipment
  - Degree of protection against harmful ingress of water: Ordinary equipment and all of applied parts (IPX7, IPX8) except ECG meet ingress protection level according to IEC 60529.
  - Degree of safety of application in the presence of a flammable anesthetic material with air or with oxygen or nitrous oxide: Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
  - Mode of operation: Continuous operation
- EN 60601-1:1988+A1:1991+A2:1995 (IEC 60601-1:1988 + A1:1991+A2:1995)
- EN 60601-1:2006+A11:2011+A1:2013 (IEC 60601-1:2005 + A1:2012)
- EN 60601-1-2:2007 (IEC 60601-1-2:2007)
- EN 60601-1-4:1996+A1:1999 (IEC 60601-1-4:1996 + A1:1999)
- EN 60601-1-6:2010 (IEC 60601-1-6:2010+A1:2013)
- EN 60601-2-37: 2008+A11:2011 (IEC 60601-2-37:2007)
- EN 55022:2010, Class B (CISPR 22:2008, modified)
- EN 55011:2009+A1:2010, Group 1, Class B
- EN 61000-4-2:2009 (IEC 61000-4-2:2008)
- EN 61000-4-3:2006+A1:2008+A2:2010 (IEC 61000-4-3:2006 + A1:2007+A2:2010)
- EN 61000-4-4:2004+A1:2010 (IEC 61000-4-4: 2004 + A1:2010)
- EN 61000-4-5:2006 (IEC 61000-4-5:2005)
- EN 61000-4-6:2009 (IEC 61000-4-6:2008)
- EN 61000-4-8:2010 (IEC 61000-4-8:2009)
- EN 61000-4-11:2004 (IEC 61000-4-11:2004)
- EN 61000-3-2:2006+A1:2009+A2:2009 (IEC 61000-3-2:2005+A1:2009+A2:2009)
- EN 61000-3-3:2008 (IEC 61000-3-3:2008)
- EN 62304:2006 (IEC 62304:2006)
- EN 62366:2008 (IEC 62366:2014)
- EN ISO 14971:2012 (ISO 14971:2007)
- EN ISO 10993-1:2009 (ISO 10993-1:2009)
- EN ISO 10993-5:2009 (ISO 10993-5:2009)
- EN ISO 10993-10:2013 (ISO 10993-10:2010)
- AIUM/NEMA UD 2:2004 (R2009)
- AIUM/NEMA UD 3:2004 (R2009)
- MDD 93/42/EEC
- CANADA [Regulation] SOR-98-282
- CE

© 2015 ALPINION MEDICAL SYSTEMS – All rights reserved.

2015 ALPINION MEDICAL SYSTEMS reserves the right to make changes in specifications and features shown hereby, or discontinue the product described at any time without notice or obligation. Contact you ALPINION representative or the most current information.

[WWW.ALPINION.COM](http://WWW.ALPINION.COM)